


SALSA® MLPA®

Cilia Disorders

Coffalyser.Net: MLPA® analysis software



Free MLPA data analysis software designed and supported by MRC Holland.

- User-friendly software and reliable MLPA data analysis
- Extensive quality control developed specifically for MLPA
- Immediate access to the latest analysis panels (Coffalyser sheets)
- Server-client model that allows data sharing
- Available free of charge!

Collaborations with scientists

Most novel MLPA applications are developed in close collaboration with scientists around the world. Results obtained with MLPA probemixes have been described in thousands of scientific publications. Researchers are encouraged to contact us with requests for new MLPA applications or feedback on current panels on info@mrcholland.com.



MLPA® & Cilia Disorders

Multiplex Ligation-dependent Probe Amplification (MLPA) is a multiplex PCR-based method that can detect the copy number of up to 60 DNA sequences in a single reaction. 96 DNA samples can be handled simultaneously, with results being available within 24 hours.

In addition to copy number changes, MLPA allows for the detection of select known point mutations. Furthermore, MLPA is able to detect methylation patterns in DNA when used in combination with a methylation-sensitive restriction enzyme (MS-MLPA). MLPA is used worldwide for diagnostics and research of human genetic disorders and tumours.

				
<p>Simultaneous detection</p> <p>of copy number, methylation and select known point mutations.</p>	<p>Low input</p> <p>Requires only 50 ng of DNA.</p>	<p>Time-efficient</p> <p>Results available within 24 hours.</p>	<p>Short hands-on time</p> <p>MLPA is performed in 5 simple steps.</p>	<p>Cost-effective</p> <p>One MLPA reaction costs EUR 12/USD 15.</p>

MLPA® Protocol

1. Sample DNA denaturation

- Sample DNA is heated to fully denature the DNA

2. Hybridisation of probes to sample DNA

- SALSA MLPA Buffer and a SALSA MLPA Probemix consisting of up to 60 probes are added to the sample

3. Ligation of hybridised probes

- Hybridised probes are ligated by adding SALSA Ligase-65 enzyme and SALSA Ligase Buffers to form fully amplifiable probes

4. PCR amplification of ligated probes






- Ligated MLPA probes are amplified by adding SALSA Polymerase and a single fluorescently-labelled primer pair

5. Fragment separation by capillary electrophoresis

- MLPA PCR products are directly loaded onto a CE device

6. Analysis by Coffalyser.Net

How MLPA® works



MLPA® Probemixes: Cilia Disorders

Over 350 MLPA probemixes are available and new assays are continuously developed in close collaboration with scientists around the world. The following lists give an overview of current MLPA probemixes for cilia disorders. See www.mrcholland.com for a complete overview.

Application	Probemix	Genes/region
Autosomal Dominant Polycystic Kidney Disease (ADPKD)	P351	<i>PKD1</i>
	P352	<i>PKD1, PKD2</i>
Autosomal Recessive Kidney Disease (ARPKD)	P341	<i>PKHD1</i>
	P342	<i>PKHD1</i>
Birt-Hogg-Dube Syndrome	P256	<i>FLCN</i>
Leber Congenital Amaurosis (LCA)	P222	<i>CEP290, GUCY2D, RDH12, RPGRIP1</i>
Nephronophthisis 1	P387	<i>NPHP1</i>
Primary Ciliary Dyskinesia (PCD) (Kartagener Syndrome, Immotile Ciliary Syndrome)	P237	<i>DNAI1</i>
	P238	<i>DNAH5</i>
Retinitis Pigmentosa (RP)	P235	<i>IMPDH1, PRPF31, RHO, RP1</i>
	P366	<i>CHM, RP2, RPGR</i>
Usher Syndrome	P292	<i>PCDH15</i>
	P361	<i>USH2A</i>
	P362	<i>USH2A</i>
Von Hippel-Lindau Syndrome	P016	<i>VHL</i>

MLPA probemixes are for Research Use Only. Not for Use in Diagnostic Procedures unless explicitly stated otherwise.

