


SALSA[®] MLPA[®]

Predisposition to Cancer

Coffalyser.Net: MLPA[®] analysis software



Free MLPA data analysis software designed and supported by MRC Holland.

- User-friendly software and reliable MLPA data analysis
- Extensive quality control developed specifically for MLPA
- Immediate access to the latest analysis panels (Coffalyser sheets)
- Server-client model that allows data sharing
- Available free of charge!

Collaborations with scientists

Most novel MLPA applications are developed in close collaboration with scientists around the world. Results obtained with MLPA probemixes have been described in thousands of scientific publications. Researchers are encouraged to contact us with requests for new MLPA applications or feedback on current panels on info@mrcholland.com.



MLPA® & Predisposition to Cancer

Multiplex Ligation-dependent Probe Amplification (MLPA) is a multiplex PCR-based method that can detect the copy number of up to 60 DNA sequences in a single reaction. 96 DNA samples can be handled simultaneously, with results being available within 24 hours.

In addition to copy number changes, MLPA allows for the detection of select known point mutations. Furthermore, MLPA is able to detect methylation patterns in DNA when used in combination with a methylation-sensitive restriction enzyme (MS-MLPA). MLPA is used worldwide for diagnostics and research of human genetic disorders and tumours.




Simultaneous detection

of copy number, methylation and select known point mutations.




Low input

Requires only 50 ng of DNA.




Time-efficient

Results available within 24 hours.



Short hands-on time

MLPA is performed in 5 simple steps.



Cost-effective

One MLPA reaction costs EUR 12/USD 15.

MLPA® Protocol

1. Sample DNA denaturation

- Sample DNA is heated to fully denature the DNA

2. Hybridisation of probes to sample DNA

- SALSA MLPA Buffer and a SALSA MLPA Probemix consisting of up to 60 probes are added to the sample

3. Ligation of hybridised probes

- Hybridised probes are ligated by adding SALSA Ligase-65 enzyme and SALSA Ligase Buffers to form fully amplifiable probes

4. PCR amplification of ligated probes

- Ligated MLPA probes are amplified by adding SALSA Polymerase and a single fluorescently-labelled primer pair


5. Fragment separation by capillary electrophoresis

- MLPA PCR products are directly loaded onto a CE device


6. Analysis by Coffalyser.Net

How MLPA® works


1. Denaturation/2. Hybridisation: Left (LPO) and Right Probe Oligo (RPO) bind to their target DNA.




3. Ligation: Hybridised probe oligos are ligated by ligase enzyme.



4. Amplification: Ligated probes are amplified using a single primer pair.




5. Fragment Separation: PCR products are separated by length.



6. Analysis and Reporting: Coffalyser.Net performs a quality check and calculates probe ratios.

A probe ratio of 1.0 signifies a normal diploid copy number; a probe ratio of 0.5 a heterozygous deletion.



MLPA® Probemixes: Predisposition to Cancer

Over 350 MLPA probemixes are available and new assays are continuously developed in close collaboration with scientists around the world. The following list gives an overview of current MLPA probemixes for predisposition to cancer. See www.mrcholland.com for a complete overview.

Application	Probemix	Main target genes			
Ataxia Telangiectasia	P041-P042	ATM			
Cowden Syndrome	P225	PTEN			
DICER1 Syndrome and Pleuro-pulmonary Blastoma	P482	DICER1			
Familial Hyperparathyroidism + HPT Jaw Tumour Syndrome	P466	CDC73			
Familial MDS-AML	P437	GATA2 (+R398W, T354M), TERC, TERT, CEBPA, RUNX1			
Familial Medulloblastoma and Meningioma	P472	SUFU			
Familial Melanoma	P419	CDKN2A/2B, CDK4, MITF E318K			
	ME024	CDKN2A/2B*, flanking regions			
Familial Meningioma	P478	SMARCE1			
Fanconi Anemia	P031-P032	FANCA			
	P057	FANCD2, PALB2			
	P113	FANCB			
	P260	PALB2, RAD50, RAD51C, RAD51D			
Gorlin Syndrome	P067	PTCH1			
Hereditary Breast & Ovarian Cancer	P002	BRCA1	Primary test	P087	Confirmatory
	P045	BRCA2, CHEK2	Primary test	P077	Confirmatory
	P090	BRCA2	Primary test	P077	Confirmatory
	P239	BRCA1 region			
Breast Cancer, Increased Susceptibility To	P041-P042	ATM			
	P057	FANCD2, PALB2			
	P190	CHEK2, ATM, TP53			
	P240	BRIP1, CHEK1			
	P260	PALB2, RAD50, RAD51D, RAD51C			
Hereditary Diffuse Gastric Cancer	P083	CDH1			
Juvenile Polyposis Syndrome	P158	BMPRIA, SMAD4, PTEN			
Li-Fraumeni Syndrome	P056	TP53, CHEK2 1100delC			

Application	Probemix	Main target genes			
Lynch Syndrome	P003	MLH1, MSH2	Primary test	P248	Confirmatory
	P008	PMS2			
	P072	MSH6, MUTYH, MSH2, EPCAM			
	ME011	EPCAM, BRAF V600E + Methylation profiling of MLH1, MSH2, MSH6, PMS2			
Melanocytic Tumours, Mesothelioma	P417	BAP1			
Multiple Endocrine Neoplasia	P017	MEN1			
	P244	AIP, MEN1, CDKN1B			
Multiple Osteochondromas	P215	EXT1, EXT2			
Neurofibromatosis	P081-P082	NF1			
	P044	NF2			
	P122	NF1-area (17q11)			
	P295	SPRED1			
Papillary Renal Carcinoma	P308	MET, PTEN, LRRK2			
Paragangliomas and Pheochromocytoma	P226	SDHB, SDHC, SDHD, SDHAF1, SDHAF2			
	P429	SDHA, MAX			
Peutz-Jeghers Syndrome	P101	STK11			
Polyposis Syndrome	P043	APC			
	P378	MUTYH, GREM1, SCG5			
Retinoblastoma	P047	RB1*			
Rhabdoid Predisposition Syndrome + Schwannomatosis	P258	SMARCB1			
	P455	LZTR1			
Tuberous Sclerosis	P046	TSC2	Primary test	P337	Confirmatory
	P124	TSC1			
Von Hippel-Lindau Syndrome	P016	VHL			
Wilms' tumour, WAGR, Denys-Drash, Frasier Syndrome	P118	WT1			

MLPA probemixes are for Research Use Only. Not for Use in Diagnostic Procedures unless explicitly stated otherwise.
* For this gene/application, both copy number and DNA methylation can be determined.