

# SALSA<sup>®</sup> MS-MLPA<sup>®</sup>

A Simple Solution for Copy Number  
Determination & Methylation Profiling





## Methylation-Specific MLPA®

Methylation-Specific (MS)-MLPA® is a variant of the MLPA technique, the gold standard in copy number determination. By combining the MLPA technique with the use of the methylation-sensitive endonuclease HhaI, it is possible to detect both DNA copy number and methylation status without bisulfite treatment. With MS-MLPA, semi-quantitative methylation profiling of multiple targets is accomplished simultaneously and with ease.

### The Principle of MS-MLPA®

MS-MLPA probes hybridise to their target sequences on the sample DNA. The MS-MLPA reaction is split in two: an undigested reaction for copy number determination, and a digested reaction for methylation profiling. In the digested reaction, the HhaI enzyme is used to digest unmethylated double-stranded probe-DNA complexes.



Complexes with methylated sample DNA are not digested by HhaI. Digested probes cannot be amplified exponentially, and therefore produce no signal, while non-digested probes generate a normal signal. Methylation is displayed as the percentage of the targeted DNA sequence that was methylated in the sample.

## MS-MLPA® Probemixes


SALSA MLPA Probemix	Application
<b>ME001 and ME002</b> Tumour suppressor mix 1 and 2	Methylation profiling for 25 tumour suppressor genes (TSG).
<b>ME011</b> Mismatch Repair Genes	Mismatch repair genes methylation and associated Lynch syndrome genome changes profiling.
<b>ME012</b> MGMT-IDH1-IDH2	Targeted glioma profiling including MGMT methylation and IDH1/2 mutation detection.
<b>ME024</b> 9p21 CDKN2A/2B region	Cell cycle regulator profiling on 9p21 region associated with multiple tumour types including melanoma.
<b>ME028</b> Prader-Willi/Angelman	Imprinting disorder profiling for Prader-Willi and Angelman syndromes.
<b>ME029</b> FMR1/AFF2	Fragile X syndrome associated promoter methylation profiling for male samples.
<b>ME030</b> BWS/RSS	Imprinting disorder profiling for Beckwith-Wiedemann and Russell-Silver syndromes.
<b>ME031</b> GNAS	Imprinting disorder profiling for Albright hereditary osteodystrophy and pseudohypoparathyroidism.
<b>ME032</b> UPD7-UPD14	Imprinting disorder profiling for UPD7, RSS, UPD14 (Temple & Kagami-Ogata syndromes).
<b>ME033</b> TNDM	Imprinting disorder profiling for Transient Neonatal Diabetes Mellitus (TNDM).
<b>ME034</b> Multi-locus Imprinting	Multi-locus imprinting disturbance profiling and to distinguish maternal and paternal triploidies.
<b>ME042</b> CIMP	CpG Island Methylator Phenotype profiling (CIMP).

# MS-MLPA® Application Highlight: Imprinting Disorders

MS-MLPA® is frequently used to detect imprinting region methylation, like the Prader-Willi/Angelman (ME028 PWS/AS) and Beckwith-Wiedemann/Russell-Silver (ME030 BWS/RSS) syndromes critical regions. These regions are known to have imprinting defects, such as those caused by uniparental disomy (UPD). MS-MLPA gives reliable results for both copy number and methylation status in the regions of interest.


## ME028 – PWS/AS Copy Number & Methylation Profiling

ME028 contains over 40 probes targeting the Prader-Willi/Angelman critical region, covering *SNRPN* and *MAGEL2* to establish DNA copy number. In addition, five dedicated methylation-sensitive probes enable methylation profiling of the region.†



## ME030 – BWS/RSS Copy Number & Methylation Profiling


ME030 contains over 40 probes that target the 11p15 Beckwith-Wiedemann/Russell-Silver critical region that are used for copy number determination. Methylation status of the region is determined by eight methylation-sensitive probes that detect sequences in the H19DMR/IC1 and KvDMR/IC2 domains.†



† Sample results must be compared to results obtained on DNA samples from unaffected individuals.

\* Images are an abridged representation of the target probe results for copy number and methylation profiling. T: target probe; H: target probe with an HhaI site; M: targeted methylation site

# SALSA® MS-MLPA® Protocol



## 1. Sample DNA denaturation

Sample DNA is heated to fully denature the DNA.

## 2. Hybridisation of probes to sample DNA

SALSA MLPA Buffer and a SALSA MS-MLPA Probemix consisting of up to 60 probes are added to the sample.

## 3. Ligation & digestion of hybridised probes

At this step, your reaction is split in two. In both tubes, hybridised probes are ligated by adding SALSA Ligase-65 enzyme and SALSA Ligase Buffers to form fully amplifiable probes. HhaI enzyme is added to one tube. This results in the digestion of non-methylated probe/allele hybrids.

## 4. PCR amplification of ligated probes

Ligated complete probes are amplified by adding SALSA Polymerase and a single fluorescently-labelled primer pair.

## 5. Fragment separation by capillary electrophoresis

MS-MLPA PCR products are loaded directly onto a capillary electrophoresis device.

## 6. Data analysis by Coffalyser.Net

## Coffalyser.Net: SALSA® MS-MLPA® Analysis



Coffalyser.Net is **free** MLPA/MS-MLPA data analysis software designed and supported by MRC Holland.

- User-friendly software for reliable MLPA and MS-MLPA data analysis.
- Extensive quality control developed specifically for MLPA and MS-MLPA data.
- Immediate access to the latest analysis panels (Coffalyser sheets).
- Server-client model that allows data sharing.