

MRC Holland Support

Support > Help Centre > digitalMLPA & Coffalyser digitalMLPA > Experimental Setup
> Sample Selection & Requirements > How many samples do I need to include in a digitalMLPA experiment with no dedicated reference samples?

How many samples do I need to include in a digitalMLPA experiment with no dedicated reference samples?

This article was retrieved from support.mrcholland.com on Tuesday, 17th June 2025.

The minimum number of samples required in a digitalMLPA™ experiment without dedicated reference samples depends on several factors, including the experimental conditions, the quality of the samples, probemix characteristics and the overall reproducibility in the experiment. Therefore, the minimum number of required samples needs to be determined experimentally under your conditions.

Background

The use of dedicated reference samples is recommended for most digitalMLPA experiments. In certain situations, e.g. when a large number of samples are used and copy number changes are expected to be (very) rare, it is possible to omit dedicated reference samples and to rely on all test samples for inter-normalisation.

The reference sample quality (RSQ) score is one of the most important readouts to determine if the number of samples was high enough for reliable normalisation. The RSQ measures the variation of probes in the sample population used for normalisation. Depending on the amount of variation between samples in your experiment, it may be necessary to include more samples to avoid an RSQ warning or error. [Read more about how the RSQ is calculated.](#)

The figure below shows an example of the effect of the number of samples in an experiment without dedicated reference samples. In this experiment, an RSQ error was triggered when eight samples were included in the analysis, an RSQ warning was (just) triggered when 24 samples were included, and no warning was triggered when 48 samples were included. The inclusion of more samples leads to a narrower distribution of [corrected MAD scores](#), and reduces the chance that an RSQ warning or error is triggered due to inherent variation.



Tags

Coffalyser digitalMLPA
digitalMLPA

Related Pages

- [How do I troubleshoot RSQ warnings or errors in my digitalMLPA experiment?](#)
- [What does the RSQ warning on my Coffalyser digitalMLPA report indicate and what is it based on?](#)

Disclaimer

The information provided in this material is correct for the majority of our products. However, for certain applications, the instructions for use may differ. In the event of conflicting information, the relevant instructions for use take precedence.