

# MRC Holland Support

Support > Help Centre > digitalMLPA & Coffalyser digitalMLPA > Experimental Setup > Sequencing > What Library Prep Kit do I select for digitalMLPA in the Illumina Local Run Manager?

---

## What Library Prep Kit do I select for digitalMLPA in the Illumina Local Run Manager?

This article was retrieved from [support.mrcholland.com](https://support.mrcholland.com) on Friday, 9th May 2025.

When sequencing SALSA® digitalMLPA™ reactions, a run sheet must be prepared in the Local Run Manager of the Illumina sequencer. Instructions on how to prepare the run sheet are given in the [digitalMLPA General Protocol](#). The Library Prep Kit that is selected has little effect on the run sheet, provided that the kit allows you to select the settings specified in the protocol, and that it does not enable adapter trimming.

If the digitalMLPA library is combined with another NGS library, you can simply select the Library Prep Kit suitable for the other library. However, adapter trimming should be *disabled* to prevent possible truncation of digitalMLPA reads due to overlap with adapter sequences in the other library. If adapter trimming is enabled by the Library Prep Kit that you selected and you cannot disable it, the FASTQ file should be re-exported with adapter trimming disabled before analysis of the digitalMLPA reads.

It is also possible to import our [digitalMLPA Library Prep Kit template](#) if you have difficulty selecting a suitable Library Prep Kit.

### Important

Carefully check all settings against the [digitalMLPA General Protocol](#) before use.

To import the digitalMLPA Library Prep Kit template:

- Download the [digitalMLPA Library Prep Kit template](#).
- Navigate to the web interface of the Illumina Local Run Manager.
- Find the Library Prep Kits in your Local Run Manager. Depending on the version the exact instructions may differ. For Local Run Manager v3.0.0, click Tools > Index & Library Prep Kits; for Local Run Manager v1.3.1, click the gear icon and Library Prep Kits.
- Click the Add Library Prep Kit button and select the downloaded template.
- You should see a message that the Library Prep Kit was successfully created.

The digitalMLPA Library Prep Kit will now be available when creating a GenerateFASTQ run sheet. This template can be used for digitalMLPA sequencing runs, but not for combined sequencing runs.

## Note

The template and instructions have been prepared for Illumina Local Run Manager v1.3.1 and v3.0.0. Instructions may be different for your version of the Local Run Manager, and the template may not work on all versions. Always check run settings against the [digitalMLPA General Protocol](#) before starting your run. [Contact us](#) if you experience any difficulties.

## Tags

digitalMLPA

## Disclaimer

---

*The information provided in this material is correct for the majority of our products. However, for certain applications, the instructions for use may differ. In the event of conflicting information, the relevant instructions for use take precedence.*