

MRC Holland Support

Support > Help Centre > MLPA & Coffalyser.Net > Coffalyser.Net Data Analysis Software > Installation & Configuration > Requirements > Can Coffalyser.Net and Coffalyser digitalMLPA be installed on macOS or Linux?

Can Coffalyser.Net and Coffalyser digitalMLPA be installed on macOS or Linux?

This article was retrieved from support.mrcholland.com on Wednesday, 24th April 2024.

It is not possible to install Coffalyser.Net or Coffalyser digitalMLPA on macOS (Mac OS X) or Linux.

Background

Both programs runs on the Microsoft .NET Framework, which means that they are dependent on the Windows operating system. However, it is possible to run Windows programs on macOS or Linux systems via a virtual machine that runs Windows. For more details, see the system requirements of [Coffalyser.Net](#) or [Coffalyser digitalMLPA](#).

Tags

Coffalyser digitalMLPA

Coffalyser.Net

Related Pages

- [What are the system requirements for Coffalyser digitalMLPA?](#)
- [What are the system requirements for Coffalyser.Net?](#)

Disclaimer

The information provided in this material is correct for the majority of our products. However, for certain applications, the instructions for use may differ. In the event of conflicting information, the relevant instructions for use take precedence.