

# MRC Holland Support

Support > Help Centre > digitalMLPA & Coffalyser digitalMLPA > Coffalyser digitalMLPA Data Analysis Software > Installation & Configuration > [How do I find which version of Coffalyser digitalMLPA I have?](#)

---

## How do I find which version of Coffalyser digitalMLPA I have?

This article was retrieved from [support.mrcholland.com](https://support.mrcholland.com) on Friday, 9th May 2025.

You can find the version of a downloaded Coffalyser digitalMLPA™ software package in the copy of the release notes (PDF) in the \_Documentation folder in the location where you unpacked Coffalyser digitalMLPA. The release notes list the version of the software package as a whole and the versions of all modules. The versions of the modules used in the analysis of an experiment are also listed in the General Details and Sample Results reports.

### Modules

Coffalyser digitalMLPA consists of several modules that are bundled in a single software package: configuration, data conversion, fragment analysis, comparative analysis, extended analysis and reporting. Each module has a version number, but the software package as a whole also has a version number that is increased every time one or more of the modules are updated.

### Note

This article is about Coffalyser digitalMLPA for digitalMLPA data analysis. For information about Coffalyser.Net for conventional MLPA data analysis, see [this article](#).

### Tags

Coffalyser digitalMLPA

---

### Disclaimer

*The information provided in this material is correct for the majority of our products. However, for certain applications, the instructions for use may differ. In the event of conflicting information, the relevant instructions for use take precedence.*