

# MRC Holland Support

Support > Help Centre > digitalMLPA & Coffalyser digitalMLPA > digitalMLPA Technique > Video: What is SALSA digitalMLPA?

---

## Video: What is SALSA digitalMLPA?

This article was retrieved from support.mrcholland.com on Sunday, 11th May 2025.

### Note

This video is about digitalMLPA. For an introductory video about conventional MLPA, see [this article](#).

SALSA® digitalMLPA™ is a technique for the copy number quantification of up to 1000 target sequences. The technique combines the conventional MLPA technology with next-generation sequencing (NGS). This video briefly describes what digitalMLPA is and what steps are involved.

For a more technical overview of the digitalMLPA technique we also recommend the video [How does SALSA digitalMLPA work?](#)

Please note that subtitles are available.

In case you're unable to view the video on YouTube you can also watch it [here](#).

### Tags

digitalMLPA

E-learning

Instruction Videos

### Related Pages

---

- [Video: How does SALSA digitalMLPA work?](#)
- [Getting started with digitalMLPA](#)
- [Video: How does MLPA work?](#)

### Disclaimer

---

*The information provided in this material is correct for the majority of our products. However, for certain applications, the instructions for use may differ. In the event of conflicting information, the relevant instructions for use take precedence.*