

MRC Holland Support

Support > Help Centre > digitalMLPA & Coffalyser digitalMLPA > Coffalyser digitalMLPA Data Analysis Software > Installation & Configuration > What are the system requirements for Coffalyser digitalMLPA?

What are the system requirements for Coffalyser digitalMLPA?

This article was retrieved from support.mrcholland.com on Friday, 26th April 2024.

Note

This article is about Coffalyser digitalMLPA for digitalMLPA data analysis. For information about Coffalyser.Net for conventional MLPA data analysis, see [this article](#). Coffalyser digitalMLPA™ does not require installation. For successful operation of the software, the following requirements should be met:

- A genuine copy of Windows 10 (32- or 64-bit) with the latest updates installed.
- An up-to-date Microsoft .NET Framework (at least version 4.8).
- An Intel i5 processor or equivalent.
- 4 GB of available RAM.
- Free disk space of at least the size of the FASTQ files that are processed.

For optimal operation, we recommend:

- An Intel i7 processor or equivalent.
- 8 GB of available RAM.
- A solid state drive (SSD) with sufficient free disk space.

The use of a USB stick to store files may greatly reduce the speed of the analysis. An internet connection is not required to run Coffalyser digitalMLPA.

Tags

Coffalyser digitalMLPA

Related Pages

- [Can Coffalyser.Net and Coffalyser digitalMLPA be installed on macOS or Linux?](#)

Related Files

- [Coffalyser digitalMLPA User Manual](#)

Disclaimer

The information provided in this material is correct for the majority of our products. However, for certain applications, the instructions for use may differ. In the event of conflicting information, the relevant instructions for use take precedence.