

# MRC Holland Support

Support > Help Centre > General Information > General Product Information > Where can I find the exact sequences recognised by SALSA (digital)MLPA probes?

---

## Where can I find the exact sequences recognised by SALSA (digital)MLPA probes?

This article was retrieved from [support.mrcholland.com](https://support.mrcholland.com) on Wednesday, 8th May 2024.

Complete target sequences of most conventional MLPA probes can be found in an Excel file titled Probe sequences under the heading Other Downloads on the product page. The product description for many probemixes also contains information about the ligation site and the 24 nt surrounding the ligation site.

Complete target sequences of digitalMLPA probes are not available, but the Probe Information File (PIF) of the probemix contains partial target locations of most probes.

The exact target sequence of some (digital)MLPA probes is not available, as this is considered proprietary information.

### Note

[Contact us](#) if you need more information to verify a single-probe copy number variation or another type of unexpected result. We can also assist you in determining if an identified SNV/mutation is expected to affect a probe's signal.

### Tags

digitalMLPA

MLPA

---

### Disclaimer

*The information provided in this material is correct for the majority of our products. However, for certain applications, the instructions for use may differ. In the event of conflicting information, the relevant instructions for use take precedence.*