

MRC Holland Support

Support > Help Centre > digitalMLPA & Coffalyser digitalMLPA > Experimental Setup > Sequencing > [Which NGS systems can be used for digitalMLPA?](#)

Which NGS systems can be used for digitalMLPA?

This article was retrieved from support.mrcholland.com on Friday, 9th June 2023.

Our digitalMLPA technique is compatible with all Illumina sequencing platforms. This includes the iSeq, MiniSeq, MiSeq, NextSeq 500/550, NextSeq 1000/2000 and NovaSeq.

The digitalMLPA technique is currently not suitable for other NGS systems, including the Ion Torrent.

Tags

digitalMLPA

Disclaimer

The information provided in this material is correct for the majority of our products. However, for certain applications, the instructions for use may differ. In the event of conflicting information, the relevant instructions for use take precedence.